

BabyFlow™ nasal CPAP – reusable and disposable components



**Optimize comfort,
efficiency and
cost-effectiveness**

Optimize comfort, efficiency and cost-effectiveness



The BabyFlow™ nCPAP system incorporates current technology to help keep premature babies comfortable. The fact that no separate nCPAP device is needed streamlines workflow for hospitals. In addition to this, the system drives down investment costs through reduced equipment complexity at the bedside. Designed specifically for the Babylog 8000 ventilator series and the Evita* family ventilators. BabyFlow helps you positively impact perinatal care.

- BabyFlow is available both as reusable and disposable product
- Anatomically shaped prong and mask design
- New reusable headbands and disposable caps
- Innovative attachment accessories
- Easy application and stable nCPAP

* NIV/NeoFlow option, latest software release required.

BabyFlow Benefits

- The customer can choose between reusable and disposable BabyFlow components
- Seamless transition to nCPAP with Babylog 8000 series ventilators/ Evita* family ventilators
- Low noise level during nCPAP application supports developmental care
- Flow-independent nCPAP control
- Minimal dead space
- No separate nCPAP device needed
- Reduces equipment complexity at the bedside
- Easy training and operation



Anatomically shaped prong and mask design

- Nasal silicon prongs in four different sizes (S, M, XM, L) and silicon NeoMasks in two different sizes (S, M) assure a comfortable fit

New reusable headbands and disposable caps

- A complete line of reusable headbands is available, consisting of seven color-coded sizes
- A complete line of disposable caps is available, consisting of six color-coded sizes
- New intermediate sizes help assure the best patient fit
- Easy access during head ultrasound and IV cannulation procedures
- Each headband has Velcro attachment points and minimal seams for enhanced patient comfort and flexibility
- Reusable headbands are made of soft microfiber material (91% microfiber and 9% elastane), making them breathable, washable, and resistant to incubator humidity and warmth
- Disposable caps are made of 95% cotton and 5% elastane





Headband sizing chart

Size	Color	Head circumference		App. weight		Order no.
		cm	inches	grams	Lbs	
XS	White	21	8.3	500 – 700	1.1 – 1.5	84 18 418
S	Yellow	22	8.7	700 – 1,000	1.5 – 2.2	84 18 419
M	Red	25	9.8	1,000 – 1,300	2.2 – 2.9	84 18 420
L	Green	27	10.6	1,300 – 1,600	2.9 – 3.5	84 18 421
XL	Orange	30	11.9	1,600 – 2,400	3.5 – 5.3	84 18 422
XXL	Blue	33	13	2,400 – 3,000	5.3 – 6.6	84 18 423
XXL plus	Black	36	14.2	3,000 >	6.6 >	84 18 540

Cap sizing chart

Size	Color	Head circumference		App. weight		Order no.
		cm	inches	grams	Lbs	
S	Lemon	22	8,7	700 – 1,000	1.5 – 2.2	84 18 534
M	Light orange	25	9.8	1,000 – 1,300	2.2 – 2.9	84 18 535
L	Light salmon	27	10.6	1,300 – 1,600	2.9 – 3.5	84 18 536
XL	Lavender	30	11.9	1,600 – 2,400	3.5 – 5.3	84 18 537
XXL	Blue	33	13	2,400 – 3,000	5.3 – 6.6	84 18 538
XXL plus	Green	36	14.2	3,000 >	6.6 >	84 18 539



Ordering information

Reusable starter package 84 18 578

Consisting of:

- 2 x BabyFlow, reusable
- 1 x Prong small (pack of 10 pcs.)
- 1 x Prong medium (pack of 10 pcs.)
- 1 x Prong medium plus (pack of 10 pcs.)
- 1 x Prong large (pack of 10 pcs.)
- 1 x Headband S
- 1 x Headband M
- 1 x Headband L
- 1 x Headband XL
- 1 x Headband XXL
- 1 x NeoMask S (pack of 10 pcs.)
- 1 x NeoMask M (pack of 10 pcs.)

Each component of the starter package is also available separately

Disposable demo package 84 18 584

Consisting of:

- 5 x BabyFlow, disposable
- 1 x Prong medium (pack of 10 pcs.)
- 1 x Prong medium plus (pack of 10 pcs.)
- 1 x NeoMask M (pack of 10 pcs.)
- 1 x Cap L (pack of 5 pcs.)
- 1 x Cap XL (pack of 5 pcs.)

Each component of the demo package is also available separately

BabyFlow, disposable (pack of 20 pcs.)	84 18 582
BabyFlow, reusable	84 18 576
Adapter for prong or NeoMask	84 18 411
Silicon hoses with connectors	84 18 414
Prongs small (pack of 10 pcs.)	84 18 415
Prongs medium (pack of 10 pcs.)	84 18 416
Prongs medium plus (pack of 10 pcs.)	84 18 531
Prong large (pack of 10 pcs.)	84 18 417
NeoMask S (pack of 10 pcs.)	84 18 491
NeoMask M (pack of 10 pcs.)	84 18 490
Fixation strap XS (pack of 10 pcs.)	84 18 471
Fixation strap S (pack of 10 pcs.)	84 18 472
Fixation strap M (pack of 10 pcs.)	84 18 473
Fixation strap L (pack of 10 pcs.)	84 18 474
Fixation strap XL (pack of 10 pcs.)	84 18 475
Fixation strap XXL (pack of 10 pcs.)	84 18 476

USA:

Draeger Medical, Inc.
 3135 Quarry Road
 Telford, PA 18969
 USA
 Tel: +1-215-721-5400
 Toll-free: +1-800-437-2437
 Fax: +1-215-723-5935
 E-mail: info@draegermed.com

Canada:

Draeger Medical Canada Inc.
 120 East Beaver Creek Road Suite 104
 Richmond Hill Ontario L4B 4V1
 CANADA
 Tel: +1-905-763-3702
 Toll-free: +1-866-343-2273
 Fax: +1-905-763-1890
 E-mail: Canada.Support@draeger.com

www.draeger.com

Manufacturer:
 Dräger Medical AG & Co. KG
 23542 Lübeck, Germany

The quality management system at
 Dräger Medical AG & Co. KG is
 certified according to ISO 13485, ISO 9001
 and Annex II of Directive 93/42/EEC
 (Medical devices).